



**Real Genomics**



**HiYield™** **Gel/PCR Fragments Extraction Kit**  
*Mini Prep DNA Extraction*  
**96-well PCR Clean Up Kit**  
*Centrifuge/Vacuum System*  
**Protocol Book**

---

■ YDF100//YDF300

■ YDF96B-2//YDF96B-4//YDF96B-10

**Ver. 3.0**



# Precautions

## I) Handling Requirements

- Do not use a kit after its expiration date has passed.
- Some reagents contain the hazardous compounds guanidine thiocyanate or guanidine hydrochloride . Do not let these reagents touch your skin, eyes, or mucous membranes. If contact does occur, wash the affected area immediately with large amounts of water. If you spill the reagents, dilute the spill with water before wiping it up.
- Do not allow reagents containing guanidine thiocyanate to mix with sodium hypochlorite solution or strong acids. This mixture can produce a highly toxic gas.

## II) Laboratory Procedures

- Handle all samples and the resulting waste as if potentially infectious, using safe laboratory procedures. As the sensitivity and titer of potential pathogens in the sample material varies, the operator has to optimize pathogen inactivation by the Lysis Buffer or take appropriate measures according to local safety regulations. RBC Bioscience does not warrant that samples treated with Lysis Buffer are completely inactivated and non-infectious. After sample processing is completed, remove and autoclave all disposable plastics, if you worked with potentially infectious sample material.
- Do not eat, drink or smoke in the laboratory work area.
- Wear protective disposable gloves, laboratory coats and eye protection when handling samples and kit reagents.
- Do not use sharp or pointed objects when working with the reagent cartridge, in order to prevent damage of the sealing foil and loss of reagent.
- Do not contaminate the reagents with bacteria, virus, or ribonuclease. Use disposable pipettes and RNase-free pipette tips only to remove aliquots from reagent bottles. Use the general precautions described in the literature.
- Wash hands thoroughly after handling samples and test reagents.

## III) Waste Handling

- Discard unused reagents and waste in accordance with country, federal, state and local regulations.

# CONTENTS

## HiYield™

### Mini Prep DNA Extraction

#### **Gel/PCR Fragments Extraction Kit** ..... 1

Cat.No. YDF100//YDF300

- **Gel Extraction Protocol** ..... 3
- **PCR Clean Up Protocol** ..... 5
- **Trouble Shooting** ..... 7

### Centrifuge/Vacuum System

#### **96-well PCR Clean Up Kit** ..... 9

Cat.No. YDF96B-2//YDF96B-4//YDF96B-10

- **Centrifuge Protocol** ..... 11
- **Vacuum Protocol** ..... 13
- **Trouble Shooting** ..... 17



Mini prep DNA Extraction

# HiYield™ Gel/PCR DNA Fragments Extraction Kit

Cat.No. YDF100/YDF300



## Kit Contents

### Cat.No. YDF100

100 mini preps / kit

DF Buffer(B049).....80ml  
Wash Buffer (concentrated)\*.....25ml  
Elution Buffer.....10ml  
2ml Collection Tube.....100pcs  
DF Cloumn.....100pcs

### Cat.No. YDF300

300 mini preps / kit

DF Buffer(B050).....200ml  
Wash Buffer (concentrated)\*\*.....40ml  
Elution Buffer.....30ml  
2ml Collection Tube.....300pcs  
DF Cloumn.....300pcs

**Additional DF Buffer (cat.No. B049/B050) can be purchased separately.**

**Sample:** 200 mg Gel / 100 µl PCR Solution  
**Yields:** Gel 70 - 80 % / PCR Recovery 80 - 90 %  
**Eective Binding Capacity:** Approx 10 µg  
**Eective Primer Removal:** < 25 bp

**Operation Time:** 20 mi ns  
**Elution Volume:** 20-50µl  
**Seq. Cut-Off:** 50bp - 10kb

\* Add 100ml ethanol (96-100%) to Wash Buffer prior to initial use.

\*\* Add 160ml ethanol (96-100%) to Wash Buffer prior to initial use.

## **Description**

*The HiYield™ Gel / PCR DNA Fragments Extraction Kit is designed to recover or concentrate DNA fragments (50bp - 10Kb) from agarose gels, PCR or other enzymatic reactions. The unique dual purpose application and high yield DNA / minicolumn make this kit exceptional value. The method uses a chaotropic salt , guanidine thiocyanate to dissolve agarose gel and denature enzymes. DNA fragments in chaotropic salt solution bind to the glass fiber matrix of the spin column. Following washing off contaminants, the purified DNA fragments are eluted by addition of low salt elution buffer or water. Salts, enzymes and unincorporated nucleotides are effectively removed from reaction mixtures without phenol extraction or alcohol precipitation .*

## **Quality Control**

*The quality of HiYield™ Gel / PCR DNA Fragments Extraction Kit is tested on a lot - to - lot basis . The efficiency of DNA recovery is tested by isolation of DNA fragments of various sizes from either aqueous solution or agarose gel . The purified DNA is checked by agarose gel analysis.*

## **Reference**

*Vogelstein , B., and Gillespie , D. ( 1979) Proc . Natl. Acad. Sci. USA 76, 615.*

## **Note**

- \* DF Buffer contains guanidine thiocyanate which is harmful and an irritant agent. During operation , always wear a lab coat , disposable gloves, and protective goggles.*

## **Gel Extraction Protocol**

### **Gel Dissociation**

1. Excise the agarose gel slice containing relevant DNA Fragments and remove extra agarose to minimize gel slice.
2. Transfer up to 300mg of the gel slice into a microcentrifuge tube (not provided).
3. Add 500  $\mu$ l of DF Buffer to the sample and mix by vortexing.
4. Incubate at 55°C for 10-15 minutes until the gel slice has been completely dissolved. During incubation, invert the tube every 2-3 mins.

**For best efficiency and more precisely at excising the agarose gel slice, RBC Gel Cutter is available. Please contact your local RBC products supplier for further information.**

### **DNA Binding**

5. Place a DF Column into a Collection Tube.
6. Apply 800 $\mu$ l of the sample mixture from previous step into the DF Column.
7. Centrifuge at max. 10,000 xg (13,000 rpm) for 30 seconds.
8. Discard the flow-through and place the DF Column back in the Collection Tube.
9. If the sample mixture is more than 800  $\mu$ l, repeat this DNA Binding Step.



## **Wash**

10. Add 600µl Wash Buffer (ethanol added) to the DF Column.
11. Centrifuge at max. 10,000 xg (13,000 rpm) for 30 seconds.
12. Discard the flow-through and place the DF Column back into the Collection Tube.

**For TAE gels, proceed to step 13. For TBE gels, repeat Wash Steps 10-12. Boric Acid is difficult to remove and can affect downstream applications, therefore double wash is recommended.**

13. Centrifuge again for 2 minutes at max 10,000 xg (13,000 rpm) to dry the column matrix.

## **DNA Elution**

14. Transfer dried DF Column into a new microcentrifuge tube (not provided).
15. Add 20-50µl Elution Buffer or water to the center of the column matrix.
16. Allow to stand for 2 minutes until Elution Buffer or water is absorbed by the matrix.
17. Centrifuge for 2 minutes at full speed to elute purified DNA.

## **PCR Clean Up Protocol**

### **Sample preparation**

1. Transfer up to 100  $\mu$ l reaction product to a microcentrifuge tube (not provided).
2. Add 5 volume of DF Buffer to 1 volume of the sample and mix by vortexing.

### **DNA Binding**


3. Place a DF Column into a Collection Tube.
4. Apply the sample mixture from previous step into the DF Column.
5. Centrifuge at max. 10,000 xg (13,000 rpm) for 30 seconds.
6. Discard the flow-through and place the DF Column back in the Collection Tube.

### **Wash**

7. Add 600  $\mu$ l of Wash Buffer (ethanol added) to the DF Column.
8. Centrifuge at max. 10,000 xg (13,000 rpm) for 30 seconds.
9. Discard the flow-through and place the DF Column back in the Collection Tube.
10. Centrifuge again for 2 minutes at max. 10,000 xg (13,000 rpm) to dry the column matrix.

### **DNA Elution**

11. Transfer dried DF Column into a new microcentrifuge tube (not provided).
12. Add 20-50  $\mu$ l Elution Buffer or water to the center of the column matrix.
13. Allow to stand for 2 minutes until Elution Buffer or water is absorbed by the matrix.
14. Centrifuge for 2 minutes at max. 10,000 xg (13,000 rpm) to elute purified DNA.



<b>Trouble Shooting</b>	
<b>Problem</b>	<b>Possible Reason</b>
<b>Low recovery of DNA fragment</b>	<b>Size of DNA fragment is more than 5 kb</b>
	<b>Ineffective DNA elution</b>
	<b>Incomplete DNA elution</b>
	<b>TAE or TBE buffer is repeatedly used or of in correct pH</b>
	<b>Overloaded column with agarose</b>
<b>Poor performance in downstream applications</b>	<b>Eluted DNA carries salt residues</b>
<b>Non - Specific DNA fragment appears in eluted DNA</b>	<b>DNA fragment is denatured and becomes single - stranded</b>
	<b>Scalpel or razor blade used to excise gel is contaminated</b>
<b>Gel slice difficult to dissolve</b>	<b>Used high percentage agarose gel (&gt; 2.5%) (not recommended)</b>
	<b>Gel slice is too big (&gt; 300 mg)</b>

**Solution**

*Use elution solution preheated to 60°C*

*DNA elution may be in effective in acidic conditions. Optimal elution pH is between 7.0 - 8.5*

*Make sure min. 30 µl is applied and Elution Buffer is applied to center of membrane . Allow time for full absorption to membrane prior to centrifugation.*

*Repeated use of TAE/TBE buffers will cause pH to increase. Use fresh TAE/TBE buffer.*

*Higher recovery is attained when smaller amounts of agarose are present. Minimize the size of the gel slice. If > 300mg, split sample and use another column.*

*Wash the column twice with Wash Buffer*

*To re - anneal ss DNA, incubate tube at 95°C and cool slowly at room temp.*

*Use a new or clean scalpel or razor blade to excise the gel*

*Add double volume of buffer for gels >2.5% .  
Incubate with mixing every 1- 2 minutes until complete dissolution .*

*Use more than one column for gel slice > 300mg*

Centrifuge/Vacuum System  
**96-well PCR Clean Up Kit**  
Cat.No. YDF96B-2//YDF96B-4//YDF96B-10



## Kit Contents

### Cat.No. YDF96B-2 Validation Kit

0.35ml DNA Binding Plate.....	2 Plates
96-well 2ml Sample Plate.....	4 Plates
96-well PCR Plate.....	2 Plates
Sealing Film.....	5 sheets
Binding Buffer.....	65ml
Wash Buffer*.....	25ml
Elution Buffer.....	30 ml

### Cat.No. YDF96B-4

0.35ml DNA Binding Plate.....	4 Plates
96-well 2ml Sample Plate.....	8 Plates
96-well PCR Plate.....	4 Plates
Sealing Film.....	10 sheets
Binding Buffer.....	125ml
Wash Buffer*.....	25ml X 2 bots
Elution Buffer.....	30 ml

### Cat.No. YDF96B-10

0.35ml DNA Binding Plate.....	10 Plates
96-well 2ml Sample Plate.....	20 Plates
96-well PCR Plate.....	10 Plates
Sealing Film.....	25 sheets
Binding Buffer.....	160ml
Wash Buffer*.....	25ml X 5 bots
Elution Buffer.....	60 ml

**Validation kit includes all components for testing centrifuge and vacuum protocol system.  
Final contents will differ depending on user preference for centrifuge or vacuum protocol.**

\* Add 4 times volume of Ethanol (96-100%) to 1 volume of Wash Buffer prior to initial use.

## **Description**

*RBC Real Genomics 96-well PCR Clean Up Kit is based on a reliable convenient silica format and is compatible with both centrifuge and vacuum protocols. It contains all the necessary plates, buffers and accessories for easy use straight from the box. Using RBC Real Genomics 96-well PCR Clean Up Kit, as many as 96 amplification reactions can be processed in parallel in 25 minutes. The high quality silica membrane DNA filter plate ensures low well to well and batch to batch variability.*

## **Principle**

*RBC Real Genomic 96-well PCR Clean Up Kit has the selective binding properties of a unique silica-gel membrane. The binding buffer is optimized for the efficient recovery of single and double stranded DNA fragments from 100bp to 10kb directly from amplification reactions, and for the quantitative removal of primers under 40 bases and incorporated nucleotides.*

*The RBC YDF96B system is ideally suited for cleaning up DNA from many applications, such as PCR, cDNA synthesis reactions, or ligation reactions.*

## **Centrifuge**

*Any Medium/High Speed Centrifuge (Beckmann/Eppendorf/Qiagen/Sigma) equipped with 2 x 96 Plate Rotor.*

## **Vacuum System**

*RBC 96 Vacuum Manifold or any standard 96/384 microplate vacuum manifold. The vacuum protocol is ideal for benchtop operation and a lower cost solution than centrifuging. However, due to differences in vacuum pressure between various manifold, vacuum systems and samples, it may be necessary to optimise this step. RBC recommend a vacuum pressure starting from 10 inches Hg up to maximum 20 inches Hg to ensure all solution passes through each well. In practice all the wells should pass through solution within 3-5 minutes, the yields will show little variation even if some wells pass the solution through at varying speeds (some wells may process very rapidly). Maintain vacuum until the solution fully passes through each well.*

## Centrifuge Protocol

### Sample Preparation

1. Add 5 times volume of Binding Buffer to the sample and mix by pipetting.
  - a. For PCR sample less than 50 $\mu$ l, add Binding Buffer directly to original PCR tube.
  - b. For PCR sample larger than 50  $\mu$ l, transfer PCR sample to a 2ml Collection Plate then add 5 times volume of Binding Buffer to each well.

### DNA Binding

2. Place a 0.35ml DNA Binding Plate in a 2 ml Sample Plate.

***This whole unit must be placed in the centrifuge, the DNA Binding Plate requires the rigid 2ml Sample Plate for each centrifugation step. The DNA Binding Plate must be sealed with sealing tape or provided Sealing Film for each centrifugation step. Sealing tape or Sealing Film may be reused if required.***

3. Transfer the PCR mixture into the 0.35ml DNA Binding Plate.

***For mixture volumes larger than the capacity of the DNA Binding Plate, repeat DNA Binding steps until all the mixture has been transferred to the membrane.***

4. Seal with provided Sealing Film.
5. Centrifuge at 3,500 rpm for 5 min.
6. Discard the flow-through and return the 0.35ml DNA Binding Plate to the 2ml Sample Plate.



## **Wash**

7. Add 300  $\mu$ l of Wash Buffer (ethanol added) to each well of the 0.35ml DNA Binding Plate.
8. Centrifuge at 3,500 rpm for 5 min.
9. Discard the flow-through and return the DNA Binding Plate to the 2ml Sample Plate.
10. Add 300  $\mu$ l of Wash Buffer (ethanol added) to each well of the 0.35ml DNA Binding Plate to wash again.
11. Centrifuge at 3,500 rpm for 5 min.
12. Centrifuge for another 10 minutes (or incubate at 60°C for 5-10 minutes) to remove ethanol residue.

## **DNA Elution**

Standard elution volume is 100  $\mu$ l. If higher DNA concentration is required, reduce the Elution Buffer (30-50  $\mu$ l) to increase DNA concentration.

13. Transfer the 0.35ml DNA Binding Plate to a clean 2 ml Sample Plate.
14. Add 50-100  $\mu$ l Elution Buffer or ddH<sub>2</sub>O (pH 8.0-8.5) to the center of each well of the 0.35ml DNA Binding Plate.
15. Allow to stand for 3 minutes until the Elution Buffer or water has been absorbed by the matrix.
16. Centrifuge at 3,500 rpm for 5 min to elute purified DNA into a clean 2 ml Sample Plate.
17. Seal with Sealing Film or parafilm sealed lid and store at 4°C. For long term storage, store at -20°C.

## Vacuum Protocol

### Sample Preparation

1. Add 5 times volume of Binding Buffer to the sample and mix by pipetting.
  - a. For PCR sample less than 50 $\mu$ l, add Binding Buffer directly to original PCR tube.
  - b. For PCR sample larger than 50  $\mu$ l, transfer PCR sample to a 2ml Sample Plate then add 5 times volume of Binding Buffer to each well.

### DNA Binding

2. Place a 0.35ml DNA Binding Plate on top of the vacuum manifold. Place a waste tray or the provided 2ml Sample/ Collection Plate underneath to collect waste.
3. Transfer the PCR mixture into the 0.35ml DNA Binding Plate.

**For mixture volumes larger than the capacity of the DNA Binding Plate, repeat DNA Binding steps until all the mixture has been transferred to the membrane.**

4. Place the 0.35ml DNA Binding Plate on top of vacuum manifold tightly.
5. Set vacuum manifold to 10-20 inches Hg for 3-5 minutes until wells have emptied.

**Refer to vacuum notes above for detailed instructions.**

6. Turn off the vacuum and apply the rest of the mixture to each well of the DNA Binding Plate.
7. Set vacuum to 10-20 inches Hg for 3-5 minutes until wells have emptied.


## **Wash**

8. Add 300  $\mu$ l of Wash Buffer (ethanol added) to each well of the 0.35ml DNA Binding Plate.
9. Set vacuum to 10-20 inches Hg for 3-5 minutes until wells have emptied.
10. Add 300  $\mu$ l of Wash Buffer (ethanol added) to each well of the 0.35ml DNA Binding Plate to wash again.
11. Set vacuum to 10-20 inches Hg for 3-5 minutes until wells have emptied.
12. Set vacuum for another 10 minutes (or incubate at 60°C for 5-10 minutes) to remove ethanol residue.

## **DNA Elution**

Standard elution volume is 100  $\mu$ l. If higher DNA concentration is required, reduce the Elution Buffer (30-50  $\mu$ l) to increase DNA concentration.

13. Place the 0.35ml DNA Binding Plate on the top of the vacuum manifold. Place the provided 96-well PCR Plate on top of the 2ml Sample Plate as support plate, then both placed inside the vacuum manifold. The 0.35ml DNA Binding plate will elute directly into the 96-well PCR Plate, which are in close proximity.
14. Add 50-100  $\mu$ l Elution Buffer or ddH<sub>2</sub>O (pH 8.0-8.5) to the center of each well of the 0.35ml DNA Binding Plate.
15. Allow to stand for 3 minutes until the Elution Buffer or water has been absorbed by the matrix.
16. Set vacuum to 10-20 inches Hg for 1-2 minutes until wells have fully emptied and all elute is in the 96-well PCR Plate.
17. Seal with Sealing Film or parafilm sealed lid and store at 4°C. For long term storage, store at -20°C.




# Multi Well Vacuum Manifold Plus

Cat.No. YVM96P

## Multi Well Vacuum Manifold Plus

The multi well plate vacuum manifold is an aluminum manifold that has been designed and optimized for the vacuum filtration of the 96-well filter plates.

- \* Comes complete with the necessary O-ring and gasket. The control block includes the vacuum pressure gauge, vacuum metering valve, vacuum release valve, and the 1/4 inch hose barb for vacuum line attachment.
- \* Included with the vacuum manifold unit is a spacer block designed to accommodate standard 350µl receiver plates. The spacer block has been optimized to reduce the space between the receiver plate and the filter plate during vacuum filtration.



## Trouble Shooting

<p><b>Low Yield from DNA Elution</b></p>	<ol style="list-style-type: none"><li>1. <b>User utilized low pH elution buffer :</b> Use provided Elution Buffer, alkaline elution buffer or ddH<sub>2</sub>O (pH &gt; 7.0).</li><li>2. <b>Incorrect addition of ethanol :</b> Add ethanol according to instructions, record on buffer bottle, after addition close buffer lid tightly.</li><li>3. <b>High pressure vacuum pump:</b> Lower pressure to 10 inches Hg to increase DNA binding time.</li><li>4. <b>Insufficient incubation time for lysis:</b> Incubate for at least 20 minutes for proteinase K complete digestion.</li><li>5. <b>Insufficient mixing of lysis buffer , proteinase K and sample:</b> Shake plate vigorously for 20 secs to ensure mixing.</li><li>6. <b>DNA binding plate and Elution Buffer at room temperature for elution step:</b> Heat both in incubator to 60°C prior to elute.</li><li>7. <b>Add additional elution buffer step for low sample amounts e.g. extra 100 µl Elution buffer .</b></li></ol>
<p><b>Eluate/Well Filter is slightly disclosed</b></p>	<ol style="list-style-type: none"><li>1. <b>Incomplete lysis of blood:</b> Shake plate vigorously for 20 secs to ensure mixing of all reagents.</li><li>2. <b>Incorrect washing procedure:</b> Check ethanol was added correctly to wash buffer and wash buffer steps. Add wash step.</li><li>3. <b>Overloaded column with sample:</b> Split samples and process in two wells.</li></ol>
<p><b>Cross Contamination</b></p>	<ol style="list-style-type: none"><li>1. <b>Well to well contamination observed:</b> Insert tip ends to bottom of wells.Fix silicon cap/ tape firmly.</li><li>2. <b>Outside contamination:</b> Operate using sterile conditions and equipment. Wear gloves at all times.</li></ol>

**Poor performance of downstream application**

1. **DNA yield too low for application: Increase amount by splitting into two wells for lysis step.**
2. **DNA concentration is too low : Concentrate DNA by precipitation or vacuum.**
3. **Contamination in eluate: Increase eluate amount to 200 µl or more and check contamination prevention procedure.**
4. **RNA interfering with enzymatic procedures: Add RNase to degrade after Step 1 of Cell Lysis.**

<b>Real Genomics</b>	<b>96-well Extraction Range</b>	<b>Method</b>
<b>YPD96B</b>	<b>96-well Plasmid Kit</b>	<b>Silica Membrane</b>
<b>YDF96B</b>	<b>96-well PCR Clean Up Kit</b>	<b>Silica Membrane</b>
<b>YVN96B</b>	<b>96-well Viral Nucleic Acid Extraction Kit</b>	<b>Silica Membrane</b>
<b>RBP96B</b>	<b>96-well Genomic DNA Extraction (Blood) Kit</b>	<b>Silica Membrane</b>

**For a complete HT solution including all necessary reagents and buffers trust the Real Genomics 96-well range. Range from 2 plates (Validation), 4 plates, 10 plates, upwards.**



***RBCBioscience***